




LE-Agarose

BioNordika Agarose Type LE

BioNordika Standard Agarose is ideal for identification and recovery of roughly >1000 bp nucleic acid fragments by electrophoresis.

- DNase and RNase free
- Sharp resolution
- Strong gels with low background
- Gelling temperature (1,5%) +36°C +/- 1,5°C
- Melting temperature (1,5%) >+90°C
- Gel strength (1%) >1,200 g/cm²



Recommended size range	Final agarose conc % w/v	
	1xTAE Buffer	1xTBE Buffer
1,000 – 23,000	0,60	0,50
800 – 10,000	0,80	0,70
400 – 8,000	1,00	0,85
300 – 7,000	1,20	1,00
200 – 4,000	1,50	1,25
100 – 3,000	2,00	1,75

Part number BN-50004
Product Standard Agarose – Type LE

Size
500 g